

Anti-Chgb Rabbit Polyclonal Antibody For Wb, Elisa - P05060

Item Number: CM0026563

Introduction

CHGB Rabbit pAb (SCG1) for WB and ELISA. Detects human and mouse secretogranin-1, a neuroendocrine secretory granule protein. Suitable for research in adrenal medulla, pheochromocytoma, and cerebrospinal fluid biomarker studies.

[Learn More](#)

Parameter	Value
Product Name	CHGB Rabbit pAb
Remarks/Alias	SCG1; CHGB
Species	Human
Gene ID (Human)	1114
Gene ID	1114
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-300 of human CHGB (NP_001810.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P05060
Protein Weight	78kDa
Shipping	Ice bag