

Anti-Chi3L2 Polyclonal Antibody For Wb, Ihc-P, Elisa - Q15782

Item Number: CM0030523



Introduction

CHI3L2 Rabbit Polyclonal Antibody validated for WB, IHC-P, and ELISA. Targets human chitinase-3-like protein 2 (YKL-39), a secreted lectin with high expression in chondrocytes, ideal for cartilage biology studies.

[Learn More](#)

Parameter	Value
Product Name	CHI3L2 Rabbit pAb
Remarks/Alias	CHIL2; YKL39; YKL-39; CHI3L2
Species	Human
Gene ID (Human)	1117
Gene ID	1117
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 27-390 of human CHI3L2 (NP_003991.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	Q15782
Protein Weight	44kDa
Shipping	Ice bag