


Anti-Chia Rabbit Polyclonal Antibody For Wb And Elisa - Q9Bzp6

Item Number: CM0030168



Introduction

Rabbit polyclonal antibody targeting human acidic mammalian chitinase (CHIA), suitable for Western blot and ELISA. Cross-reacts with human, mouse, and rat. Use for research on chitin metabolism, Th2 immunity, and lung inflammation.

[Learn More](#)

Parameter	Specification
Product Name	CHIA Rabbit pAb
Remarks/Alias	CHIT2; AMCASE; TSA1902; CHIA
Species	Human
GeneID (human)	27159
GeneID	27159
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-368 of human CHIA (NP_001244932.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BZP6
Protein Weight	52kDa
Shipping	Ice bag