


Anti-Chia Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q9Bzp6

Item Number: CM0028254



Introduction

High-quality rabbit polyclonal antibody against human acidic mammalian chitinase (CHIA) for Western blot, IHC, and ELISA applications. Suitable for chitin metabolism, Th2 immunity, and asthma research. Cross-reacts with mouse and rat.

[Learn More](#)

Parameter	Value
Product Name	CHIA Rabbit pAb
Remarks/Alias	CHIT2; AMCASE; TSA1902; CHIA
Species	Human
GeneID (human)	27159
GeneID	27159
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-368 of human CHIA (NP_001244932.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9BZP6
Protein Weight	52kDa
Shipping	Ice bag