

Anti-Chicken IL12B Rabbit Polyclonal Antibody For Wb And Elisa - Q6X0K9

Item Number: CM0023560

Introduction

Chicken IL12B rabbit polyclonal antibody validated for WB and ELISA. Recognizes chicken IL-12B (IL-12p40), a 36 kDa secreted cytokine subunit. Ideal for avian immunology research.

[Learn More](#)

Parameter	Value
Product Name	Chicken IL12B Rabbit pAb
Remarks/Alias	IL-12B; IL-12p40
Species	Gallus
GeneID	404671
Immunogen	Recombinant protein corresponding to a sequence within amino acids 21-315 of Chicken IL12B (NP_998736.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Gallus
SWISS	Q6X0K9
Protein Weight	36kDa
Shipping	Ice bag