

Anti-Chit1 Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q13231

Item Number: CM0019591



Introduction

Rabbit polyclonal antibody against human CHIT1 (Chitotriosidase-1), suitable for WB, IHC-P, and ELISA. Cross-reacts with human and mouse. Ideal for chitin metabolism and immune defense research.

[Learn More](#)

Parameter	Value
Product Name	CHIT1 Rabbit pAb
Remarks/Alias	CHI3; CHIT; CHITD; CHIT1
Species	Human
GeneID (Human)	1118
GeneID	1118
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 217-466 of human CHIT1 (NP_003456.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q13231
Protein Weight	52kDa
Shipping	Ice bag