


Anti-Chk1 Polyclonal Antibody For Wb, Elisa - O14757

Item Number: CM0022767



Introduction

Rabbit polyclonal antibody against human Chk1 (Checkpoint kinase 1), validated for WB and ELISA. Suitable for cell cycle checkpoint and DNA damage response research. Recognizes human Chk1 protein, ~54 kDa.

[Learn More](#)

Parameter	Value
Product Name	Chk1 Rabbit pAb
Remarks/Alias	CHK1; Chk1
Species	Human
GeneID (Human)	1111
GeneID	1111
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 148-307 of Chk1 (NP_001265.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O14757
Protein Weight	54kDa
Shipping	Ice bag