


Anti-Chk1 Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - O14757

Item Number: CM0009279



Introduction

High-quality anti-Chk1 rabbit monoclonal antibody for Western blot, immunohistochemistry (paraffin), and ELISA. Reacts with human Chk1, a 54kDa checkpoint kinase involved in DNA damage response and cell cycle regulation. Ideal for oncology and genomic integrity studies.

[Learn More](#)

Parameter	Value
Product Name	Chk1 Rabbit mAb
Remarks/Alias	CHK1; CHEK1; Chk1
Species	Human
Gene ID (Human)	1111
Gene ID	1111
Immunogen	Synthetic Peptide A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Chk1 (O14757).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human
SWISS	O14757
Protein Weight	54kDa
Shipping	Ice bag