


Anti-Chk1 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - O14757

Item Number: CM0028010



Introduction

Rabbit anti-Chk1 polyclonal antibody validated in WB, IHC-P, IF/ICC, and ELISA. Detects human, mouse, and rat Chk1 (54 kDa). Ideal for cell cycle checkpoint and DNA damage repair research.

[Learn More](#)

Parameter	Value
Product Name	Chk1 Rabbit pAb
Remarks/Alias	CHK1; Chk1
Species	Human
Gene ID (Human)	1111
Gene ID	1111
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 148-307 of Chk1 (NP_001265.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	O14757
Protein Weight	54kDa
Shipping	Ice bag