

Anti-Chk2 (Chk2) Rabbit Monoclonal Antibody For Wb, If/Icc, Ihc-P, Elisa - O96017

Item Number: CM0003957



Introduction

Anti-Chk2 rabbit monoclonal antibody validated for WB, IF/ICC, IHC-P, ELISA. Recognizes human Chk2 (O96017), a DNA damage checkpoint kinase. Ideal for cell cycle, cancer, and apoptosis research.

[Learn More](#)

Parameter	Value
Product Name	Chk2 Rabbit mAb
Remarks/Alias	CDS1; CHK2; LFS2; RAD53; hCds1; HuCds1; PP1425; Chk2
Species	Human
Gene ID	11200
Immunogen	Synthetic Peptide · A synthetic peptide corresponding to a sequence within amino acids 1-100 of human Chk2 (O96017).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IF/ICC, IHC-P, ELISA
Cross Reactivity	Human
SWISS	O96017
Protein Weight	15-38 kDa/50-65 kDa
Shipping	Ice bag