

Anti-Chm Rabbit Monoclonal Antibody For Wb, Elisa - P24386

Item Number: CM0009168

Introduction

Anti-CHM Rabbit monoclonal antibody for Western blot and ELISA. Detects human Rab escort protein 1 (REP-1), essential for Rab protein prenylation. Suitable for cytosolic protein analysis. - P24386

[Learn More](#)

Parameter	Value
Product Name	CHM Rabbit mAb
Remarks/Alias	TCD; GGTA; REP-1; DXS540; HSD-32; CHM
Species	Human
Gene ID (Human)	1121
Gene ID	1121
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CHM (P24386)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	P24386
Protein Weight	73kDa
Shipping	Ice bag