

## Anti-Chm Rabbit Polyclonal Antibody For Wb, Elisa - P24386

Item Number: CM0026931

### Introduction

High-quality rabbit polyclonal antibody targeting human CHM (REP-1, choroideremia protein).  
Recombinant immunogen (aa 120-210).  
Validated for Western blotting and ELISA.  
Essential for Rab geranylgeranylation and trafficking research.

[Learn More](#)

| Parameter        | Value                                                                                                            |
|------------------|------------------------------------------------------------------------------------------------------------------|
| Product Name     | CHM Rabbit pAb                                                                                                   |
| Remarks/Alias    | TCD; GGTA; REP-1; DXS540; HSD-32; CHM                                                                            |
| Species          | Human                                                                                                            |
| GeneID (Human)   | 1121                                                                                                             |
| GeneID           | 1121                                                                                                             |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 120-210 of human CHM (NP_000381.1) |
| Source           | Rabbit                                                                                                           |
| Category         | Polyclonal Antibodies                                                                                            |
| Application      | WB, ELISA                                                                                                        |
| Cross Reactivity | Human                                                                                                            |
| SWISS            | P24386                                                                                                           |
| Protein Weight   | 73kDa                                                                                                            |
| Shipping         | Ice bag                                                                                                          |