

Anti-Chml Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P26374

Item Number: CM0024495

Introduction

Rabbit polyclonal antibody targeting human CHML (REP2), component A of Rab geranylgeranyltransferase, validated for WB, IF/ICC, ELISA with mouse and rat cross-reactivity. Suitable for prenylation studies.

[Learn More](#)

Parameter	Value
Product Name	CHML Rabbit pAb
Remarks/Alias	REP2; CHML
Species	Human
Human GeneID	1122
GeneID	1122
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human CHML (NP_001812.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P26374
Protein Weight	74kDa
Shipping	Ice bag