

Anti-Chmp1A Polyclonal Antibody For Wb And Elisa - Q9Hd42

Item Number: CM0016060



Introduction

High-quality rabbit polyclonal antibody against CHMP1A, validated for Western blot and ELISA. Cross-reacts with human and mouse. Ideal for studies of ESCRT-III complex, multivesicular body formation, cytokinesis, and chromatin regulation.

[Learn More](#)

Parameter	Value
Product Name	CHMP1A Rabbit pAb
Remarks/Alias	PCH8; CHMP1; PRSM1; PCOLN3; VPS46A; VPS46-1; CHMP1A
Species	Human
GeneID (Human)	5119
GeneID	5119
Immunogen	Recombinant protein corresponding to amino acids 1-196 of human CHMP1A (NP_002759.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9HD42
Protein Weight	22kDa
Shipping	Ice bag