


Anti-Chmp1B Rabbit Polyclonal Antibody For Wb, Elisa - Q7Lbr1

Item Number: CM0025625



Introduction

CHMP1B Rabbit Polyclonal Antibody for Western blot (WB) and ELISA applications. Detects human, mouse, and rat CHMP1B, a key ESCRT-III component in endosomal sorting and cytokinesis.

[Learn More](#)

Parameter	Value
Product Name	CHMP1B Rabbit pAb
Remarks/Alias	Vps46B; C10orf2; C18orf2; CHMP1.5; Vps46-2; C18-ORF2; hVps46-2; CHMP1B
Species	Human
Gene ID (Human)	57132
Gene ID	57132
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-199 of human CHMP1B (NP_065145.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q7LBR1
Protein Weight	22kDa
Shipping	Ice bag