


Anti-Chmp2A Polyclonal Antibody For Wb, Elisa - 043633

Item Number: CM0015028



Introduction

Rabbit polyclonal antibody targeting human CHMP2A, a core ESCRT-III component. Suitable for WB and ELISA applications. Ideal for studies on endosomal sorting, cytokinesis, and viral budding.

[Learn More](#)

Parameter	Value
Product Name	CHMP2A Rabbit pAb
Remarks/Alias	BC2; BC-2; VPS2; CHMP2; VPS2A; CHMP2A
Species	Human
Gene ID (Human)	27243
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-222 of human CHMP2A (NP_055268.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	O43633
Protein Weight	25kDa
Shipping	Ice bag