

Anti-Chmp2B Polyclonal Antibody For Wb, If/Icc, Elisa - Q9Uqn3

Item Number: CM0014062



Introduction

Rabbit polyclonal antibody against human CHMP2B, suitable for WB, IF/ICC, ELISA. Cross-reacts with human and mouse. Target protein is a core component of ESCRT-III involved in MVB sorting and membrane fission.

[Learn More](#)

Parameter	Value
Product Name	CHMP2B Rabbit pAb
Remarks/Alias	DMT1; ALS17; VPS2B; VPS2-2; CHMP2.5; FTDALS7; CHMP2B
Species	Human
Gene ID (Human)	25978
Gene ID	25978
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-213 of human CHMP2B (NP_054762.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q9UQN3
Protein Weight	24kDa
Shipping	Ice bag