

Anti-Chmp2B Rabbit Monoclonal Antibody For Wb, Ihc-P, Elisa - Q9Uqn3

Item Number: CM0003233



Introduction

CHMP2B rabbit monoclonal antibody validated for Western blot, immunohistochemistry, and ELISA. Suitable for human, mouse, rat samples. Target is ESCRT-III complex protein involved in endosomal sorting and neurodegeneration.

[Learn More](#)

Parameter	Value
Product Name	CHMP2B Rabbit mAb
Remarks/Alias	DMT1; ALS17; VPS2B; VPS2-2; CHMP2.5; FTDALS7; CHMP2B
Species	Human
GeneID (human)	25978
GeneID	25978
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 114-213 of human CHMP2B (Q9UQN3).
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UQN3
Protein Weight	24kDa
Shipping	Ice bag