

Anti-Chmp4B Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q9H444

Item Number: CM0027374



Introduction

Rabbit polyclonal antibody targeting CHMP4B (VPS32B), validated for WB, IF/ICC, and ELISA. Reacts with human, mouse, and rat. Essential for ESCRT-III complex studies in MVB formation, cytokinesis, and virology research.

[Learn More](#)

Parameter	Description
Product Name	CHMP4B Rabbit pAb
Remarks/Alias	SNF7; CTPP3; Shax1; CHMP4A; SNF7-2; VPS32B; CTRCT31; Vps32-2; C20orf178; dj553F4.4; CHMP4B
Species	Human
GeneID (human)	128866
GeneID	128866
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9H444
Protein Weight	25kDa
Shipping	Ice bag