

Anti-Chn1 Rabbit Polyclonal Antibody For Wb And Elisa - P15882

Item Number: CM0025682



Introduction

CHN1 Rabbit Polyclonal Antibody validated for Western blot and ELISA, reacting with Human, Mouse, and Rat. Raised against recombinant CHN1 protein, ideal for investigating neuronal development and axon guidance.

[Learn More](#)

Parameter	Value
Product Name	CHN1 Rabbit pAb
Remarks/Alias	NC, CHN, DURS2, ARHGAP2, RHOGAP2, CHN1
Species	Human
GeneID	1123
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-200 of human CHN1 (NP_001813.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P15882
Protein Weight	53kDa
Shipping	Ice bag