

Anti-Chodl Rabbit Polyclonal Antibody For Wb, Elisa - Q9H9P2

Item Number: CM0020876

Introduction

Rabbit polyclonal antibody targeting human CHODL (chondrolectin) for WB and ELISA applications. Detects mouse and rat isoforms. Suitable for research on neurite outgrowth, axon guidance, and nervous system development.

[Learn More](#)

Parameter	Value
Product Name	CHODL Rabbit pAb
Remarks/Alias	MT75; PRED12; C21orf68; CHODL
Species	Human
GeneID (Human)	140578
GeneID	140578
Immunogen	Recombinant protein Recombinant fusion protein containing a sequence corresponding to amino acids 110-210 of human CHODL (NP_001191103.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9H9P2
Protein Weight	30kDa
Shipping	Ice bag