


Anti-Chpt1 Polyclonal Antibody For Wb, Elisa - Q8Wud6

Item Number: CM0021104



Introduction

Rabbit polyclonal antibody against human CHPT1, validated for WB and ELISA.

Cholinephosphotransferase 1 catalyzes the final step of phosphatidylcholine synthesis, crucial for membrane biogenesis.

[Learn More](#)

Parameter	Value
Product Name	CHPT1 Rabbit pAb
Remarks/Alias	CPT, CPT1, CHPT1
Species	Human
Gene ID	56994
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 1-100 of human CHPT1 (NP_064629.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q8WUD6
Protein Weight	45kDa
Shipping	Ice bag