

Anti-Chrac1 Polyclonal Antibody For Wb, Elisa - Q9Nrg0

Item Number: CM0013408

Introduction

Rabbit polyclonal antibody against chromatin accessibility complex protein 1 (CHRAC1). Validated for Western blot and ELISA, recognizes human, mouse, and rat CHRAC1. Suitable for chromatin remodeling and epigenetics studies.

[Learn More](#)

Parameter	Value
Product Name	[KO Validated] CHRAC1 Rabbit pAb
Remarks/Alias	YCL1; CHARC1; CHARC15; CHRAC-1; CHRAC15; CHRAC-15; C1
Species	Human
Gene ID (Human)	54108
Gene ID	54108
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-131 of human CHRAC1 (NP_059140.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9NRG0
Protein Weight	15kDa
Shipping	Ice bag