

Anti-Chrac1 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Nrg0

Item Number: CM0029258

Introduction

CHRAC1 Rabbit Polyclonal Antibody for WB and ELISA, targeting human chromatin accessibility complex protein 1 (CHRAC1, 15 kDa). Validated for research use in chromatin remodeling and DNA replication studies.

[Learn More](#)

Parameter	Specification
Product Name	CHRAC1 Rabbit pAb
Remarks/Alias	YCL1; CHARC1; CHARC15; CHRAC-1; CHRAC15; CHRAC-15; CHRAC1
Species	Human
GeneID (Human)	54108
GeneID	54108
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-131 of human CHRAC1 (NP_059140.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9NRG0
Protein Weight	15kDa
Shipping	Ice bag