

Anti-Chrd Rabbit Polyclonal Antibody For Wb, Elisa - Q9H2X0

Item Number: CM0022149

Introduction

Rabbit polyclonal antibody against human CHRD (Chordin), validated for WB and ELISA. Suitable for cross-reactivity with mouse and rat. Ideal for developmental biology and BMP signaling research.

[Learn More](#)

Parameter	Value
Product Name	CHRD Rabbit pAb
Remarks/Alias	CHRD
Species	Human
GeneID (Human)	8646
GeneID	8646
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 686-955 of human CHRD (NP_003732.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	Q9H2X0
Protein Weight	102kDa
Shipping	Ice bag