

Anti-Chrfam7A Rabbit Polyclonal Antibody For Wb, Elisa - Q494W8

Item Number: CM0014776

Introduction

Rabbit polyclonal antibody targeting human CHRFAM7A fusion protein. Validated for Western blot and ELISA, cross-reactive with mouse and rat. Ideal for neuroscience studies of hippocampal membrane protein expression.

[Learn More](#)

Parameter	Value
Product Name	CHRFAM7A Rabbit pAb
Remarks/Alias	D-10; CHRNA7; NACHRA7; CHRNA7-DR1; CHRFAM7A
Species	Human
Gene ID	89832
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human CHRFAM7A (NP_647536.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q494W8
Protein Weight	46kDa
Shipping	Ice bag