

Anti-Chrfam7A Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - Q494W8

Item Number: CM0028805

Introduction

Rabbit polyclonal antibody against human CHRFA7A (CHRNA7-FAM7A fusion protein), validated for Western blot, IHC-P, and ELISA. Cross-reacts with mouse and rat. Ideal for neuroscience studies.

[Learn More](#)

Parameter	Specification
Product Name	CHRFAM7A Rabbit pAb
Remarks/Alias	D-10; CHRNA7; NACHRA7; CHRNA7-DR1; CHRFAM7A
Species	Human
GeneID (Human)	89832
GeneID	89832
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-150 of human CHRFAM7A (NP_647536.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q494W8
Protein Weight	46kDa
Shipping	Ice bag