

Anti-Chrm1 Polyclonal Antibody For Wb, If/Icc, Elisa - P11229

Item Number: CM0022849



Introduction

CHRM1 Rabbit Polyclonal Antibody targeting muscarinic acetylcholine receptor M1 (P11229). Validated in WB, IF/ICC, ELISA across human, mouse, rat. Ideal for GPCR and neuroscience studies.

[Learn More](#)

Parameter	Description
Product Name	CHRM1 Rabbit pAb
Remarks/Alias	M1, HM1, M1R, CHRM1
Species	Human
GeneID (Human)	1128
GeneID	1128
Immunogen	Synthetic peptide
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11229
Protein Weight	51kDa
Shipping	Ice bag