

Anti-Chrm1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P11229

Item Number: CM0021884



Introduction

CHRM1 Rabbit pAb is a polyclonal antibody targeting muscarinic acetylcholine receptor M1, validated for WB, IHC-P, IF/ICC, and ELISA in human, mouse, and rat samples. Ideal for GPCR signaling and neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	CHRM1 Rabbit pAb
Remarks/Alias	M1, HM1, M1R, CHRM1
Species	Human
Gene ID (Human)	1128
Gene ID	1128
Immunogen	Synthetic peptide A synthetic peptide corresponding to a sequence within amino acids 251-350 of human CHRM1 (NP_000729.2).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P11229
Protein Weight	51kDa
Shipping	Ice bag