


## Anti-Chrm3 Polyclonal Antibody For Wb, Elisa - P20309

Item Number: CM0022004



### Introduction

Rabbit polyclonal antibody against human CHRM3. Suitable for WB and ELISA, cross-reacts with human, mouse, rat. Recombinant immunogen (aa 390-475). Ideal for GPCR signaling studies.

[Learn More](#)

| Parameter        | Value                                                                                                              |
|------------------|--------------------------------------------------------------------------------------------------------------------|
| Product Name     | CHRM3 Rabbit pAb                                                                                                   |
| Remarks/Alias    | HM3, PBS, EGBRS, CHRM3                                                                                             |
| Species          | Human                                                                                                              |
| GeneID (Human)   | 1131                                                                                                               |
| GeneID           | 1131                                                                                                               |
| Immunogen        | Recombinant fusion protein containing a sequence corresponding to amino acids 390-475 of human CHRM3 (NP_000731.1) |
| Source           | Rabbit                                                                                                             |
| Category         | Polyclonal Antibodies                                                                                              |
| Application      | WB, ELISA                                                                                                          |
| Cross Reactivity | Human, Mouse, Rat                                                                                                  |
| SWISS            | P20309                                                                                                             |
| Protein Weight   | 66kDa                                                                                                              |
| Shipping         | Ice bag                                                                                                            |