


Anti-Chrm5 Polyclonal Antibody For Wb And Elisa - P08912

Item Number: CM0024980



Introduction

CHRM5 Rabbit pAb validated for WB and ELISA, reacting with human, mouse, and rat. Ideal for GPCR signaling research and muscarinic receptor studies in neuroscience and pharmacology.

[Learn More](#)

Parameter	Value
Product Name	CHRM5 Rabbit pAb
Remarks/Alias	HM5; CHRM5
Species	Human
GeneID (human)	1133
GeneID	1133
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 275-370 of human CHRM5 (NP_036257.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P08912
Protein Weight	60kDa
Shipping	Ice bag