

Anti-Chrna1 Polyclonal Antibody For Wb, If, Ihc, Elisa - P02708

Item Number: CM0026988



Introduction

High-quality rabbit polyclonal antibody against human CHRNA1, suitable for WB, IF-P, IHC-P, and ELISA. Detects human and mouse CHRNA1 protein, ~52kDa. Ideal for neuromuscular and ion channel research.

[Learn More](#)

Parameter	Value
Product Name	CHRNA1 Rabbit pAb
Remarks/Alias	ACHRA; ACHRD; CHRNA; CMS1A; CMS1B; CMS2A; FCCMS; SCCMS; CHRNA1
Species	Human
GeneID (Human)	1134
Immunogen	Recombinant protein; Recombinant fusion protein containing a sequence corresponding to amino acids 105-220 of human CHRNA1 (NP_000070.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	P02708
Protein Weight	52kDa
Shipping	Ice bag