


Anti-Chrna3 Rabbit Polyclonal Antibody For Wb And Elisa - P32297

Item Number: CM0024432



Introduction

High-quality CHRNA3 rabbit polyclonal antibody validated for Western blot and ELISA applications. Cross-reacts with human, mouse, and rat samples. Ideal for neuroscience and receptor research.

[Learn More](#)

Parameter	Value
Product Name	CHRNA3 Rabbit pAb
Remarks/Alias	LNCR2; PAOD2; BAIPRCK; NACHRA3; CHRNA3
Species	Human
Gene ID (Human)	1136
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 32-240 of human CHRNA3 (NP_000734.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P32297
Protein Weight	57kDa
Shipping	Ice bag