

Anti-Chrna5 Rabbit Polyclonal Antibody For Wb And Elisa - P30532

Item Number: CM0021190

Introduction

Rabbit polyclonal antibody targeting CHRNA5, validated for Western blotting and ELISA. Reacts with human and mouse CHRNA5, a neuronal nicotinic receptor subunit involved in nicotine addiction. Suitable for neuroscience and addiction research.

[Learn More](#)

Parameter	Value
Product Name	CHRNA5 Rabbit pAb
Remarks/Alias	LNCR2; CHRNA5
Species	Human
Gene ID (Human)	1138
Gene ID	1138
Immunogen	Recombinant fusion protein containing amino acids 23-240 of human CHRNA5 (NP_000736.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	P30532
Protein Weight	53kDa
Shipping	Ice bag