

Anti-Chrna6 Rabbit Polyclonal Antibody For Wb, Elisa - Q15825

Item Number: CM0028594



Introduction

High-quality rabbit polyclonal antibody against human CHRNA6, validated for WB and ELISA, with cross-reactivity to mouse. Ideal for studying nicotinic acetylcholine receptors in neuroscience research.

[Learn More](#)

Parameter	Value
Product Name	CHRNA6 Rabbit pAb
Remarks/Alias	CHNRA6; CHRNA6
Species	Human
GeneID (human)	8973
GeneID	8973
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q15825
Protein Weight	57kDa
Shipping	Ice bag