

Anti-Chrna7 Rabbit Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - P36544

Item Number: CM0022961



Introduction

High-quality rabbit polyclonal antibody against CHRNA7 (alpha-7 nAChR). Validated for WB, IHC-P, IF/ICC, ELISA. Reacts with human, mouse, rat. Recombinant immunogen targets extracellular domain. Ideal for neuroinflammation and cholinergic signaling research.

[Learn More](#)

Parameter	Value
Product Name	CHRNA7 Rabbit pAb
Remarks/Alias	NACHRA7; CHRNA7-2; CHRNA7
Species	Human
Gene ID (Human)	1139
Gene ID	1139
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 23-230 of human CHRNA7 (NP_001177384.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P36544
Protein Weight	54kDa
Shipping	Ice bag