


Anti-Chrna9 Rabbit Polyclonal Antibody For Wb, Elisa - Q9Ugm1

Item Number: CM0025482



Introduction

High-quality rabbit polyclonal antibody targeting human neuronal acetylcholine receptor subunit alpha-9 (CHRNA9). Validated for WB and ELISA, with reactivity across human, mouse, and rat samples. Ideal for neuroscience and ion channel research.

[Learn More](#)

Parameter	Value
Product Name	CHRNA9 Rabbit pAb
Remarks/Alias	NACHRA9; HSA243342; CHRNA9
Species	Human
GeneID (Human)	55584
GeneID	55584
Immunogen	Recombinant protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q9UGM1
Protein Weight	55kDa
Shipping	Ice bag