


Anti-Chrnb2 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - P17787

Item Number: CM0030557



Introduction

Rabbit polyclonal antibody targeting human CHRNB2 (neuronal nicotinic receptor subunit beta-2). Suitable for WB, IF/ICC, ELISA with cross-reactivity to mouse and rat. Ideal for neuroscience and ion channel research.

[Learn More](#)

Parameter	Value
Product Name	CHRN2 Rabbit pAb
Remarks/Alias	EFNL3; nAChRB2; CHRN2
Species	Human
GeneID (Human)	1141
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 340-430 of human CHRN2 (NP_000739.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P17787
Protein Weight	57kDa
Shipping	Ice bag