

Anti-Chrne Rabbit Polyclonal Antibody For Wb, Elisa - Q04844

Item Number: CM0016215



Introduction

CHRNE rabbit polyclonal antibody validated for WB and ELISA, recognizing human, mouse, rat acetylcholine receptor subunit epsilon. Ideal for neuromuscular junction and ion channel research. SWISS: Q04844.

[Learn More](#)

Parameter	Value
Product Name	CHRNE Rabbit pAb
Remarks/Alias	MGI; FIM1; FIMG; ACHRE; CMS1D; CMS1E; CMS2A; CMS4A; CMS4B; CMS4C; FCCMS; FIMG1; SCCMS; CMS1A1; CHRNE
Species	Human
GeneID (Human)	1145
GeneID	1145
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 21-239 of human CHRNE (NP_000071.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q04844
Protein Weight	55kDa
Shipping	Ice bag