

Anti-Chromogranin A Rabbit Monoclonal Antibody For If-P, Ihc-P, Elisa - P10645

Item Number: CM0010910



Introduction

Chromogranin A Rabbit Monoclonal Antibody validated for IF-P, IHC-P, ELISA. Cross-reacts with human, mouse, rat. Ideal for neuroendocrine, cardiovascular, and innate immunity research. - P10645

[Learn More](#)

Parameter	Value
Product Name	Chromogranin A Rabbit mAb
Remarks/Alias	CGA; PHE5; PHES; Chromogranin A
Species	Human
GeneID	1113
Immunogen	Synthetic peptide corresponding to a sequence within amino acids 1-100 of human Chromogranin A (P10645)
Source	Rabbit
Category	Monoclonal Antibodies
Application	IF-P, IHC-P, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	P10645
Protein Weight	51kDa
Shipping	Ice bag