

Anti-Chromogranin B Rabbit Monoclonal Antibody For Wb, Elisa - P05060

Item Number: CM0010923

Introduction

Chromogranin B (SCG1) rabbit monoclonal antibody for WB and ELISA, cross-reacting with human and rat. Ideal for neuroendocrine secretory granule research and diagnostic development.

[Learn More](#)

Parameter	Value
Product Name	Chromogranin B Rabbit mAb
Remarks/Alias	SCG1; Chromogranin B
Species	Human
GenelD (Human)	1114
GenelD	1114
Immunogen	Synthetic Peptide
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human, Rat
SWISS	P05060
Protein Weight	78 kDa
Shipping	Ice bag