

Anti-Chtf18 Rabbit Polyclonal Antibody For Wb, If/Icc, Elisa - Q8Wvb6

Item Number: CM0024489



Introduction

Rabbit polyclonal antibody against human CHTF18, validated for WB, IF/ICC, ELISA. Recognizes human and rat CHTF18, a chromosome cohesion factor. Ideal for cell cycle and DNA replication studies.

[Learn More](#)

Parameter	Value
Product Name	CHTF18 Rabbit pAb
Remarks/Alias	CHL12; Ctf18; RUVBL; C16orf41; C321D2.2; C321D2.3; C321D2.4; CHTF18
Species	Human
GeneID (human)	63922
GeneID	63922
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 794-975 of human CHTF18 (NP_071375.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, ELISA
Cross Reactivity	Human, Rat
SWISS	Q8WVB6
Protein Weight	107kDa
Shipping	Ice bag