

Anti-Chuk (Ikka) Polyclonal Antibody For If/Icc, If-P, Elisa - G3V926

Item Number: CM0020569

Introduction

Rabbit polyclonal antibody targeting rat CHUK (IKK α , I κ BKA). Validated for IF/ICC, IF-P, ELISA with cross-reactivity for human and mouse. Ideal for NF- κ B signaling research.

[Learn More](#)

Parameter	Value
Product Name	CHUK Rabbit pAb
Remarks/Alias	Ikka; Ikbka
Species	Rat
GeneID (Human)	1147
GeneID	309361
Immunogen	Recombinant protein corresponding to a sequence within amino acids 400-660 of rat CHUK(NP_001101058.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	IF/ICC, IF-P, ELISA
Cross Reactivity	Human, Mouse
SWISS	G3V926
Protein Weight	78kDa
Shipping	Ice bag