

Anti-Chx10 (Vsx2) Rabbit Monoclonal Antibody For Wb And Elisa - P58304

Item Number: CM0002269

Introduction

High-specificity rabbit monoclonal antibody targeting human CHX10 (VSX2), validated for Western blot and ELISA with cross-reactivity to mouse and rat. Ideal for retinal development and transcription research.

[Learn More](#)

Parameter	Value
Product Name	CHX10 Rabbit mAb
Remarks/Alias	RET1; CHX10; HOX10; MCOP2; MCOPCB3
Species	Human
GeneID (Human)	338917
GeneID	338917
Immunogen	A synthetic peptide corresponding to a sequence within amino acids 262-361 of human CHX10 (P58304)
Source	Rabbit
Category	Monoclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse, Rat
SWISS	P58304
Protein Weight	39kDa
Shipping	Ice bag