


Anti-Ciap1 Polyclonal Antibody For Wb, Ihc-P, If/Icc, Elisa - Q13490

Item Number: CM0017475



Introduction

Rabbit polyclonal antibody detecting human, mouse, and rat cIAP1 (BIRC2). Validated for Western blot, IHC-P, IF/ICC, and ELISA. Targets E3 ubiquitin ligase regulating apoptosis and NF-κB signaling.

[Learn More](#)

Parameter	Value
Product Name	cIAP1 Rabbit pAb
Remarks/Alias	API1, MIHB, HIAP2, RNF48, cIAP1, Hiap-2, c-IAP1
Species	Human
GeneID (Human)	329
GeneID	329
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 331-420 of human cIAP1 (NP_001157.1).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, IF/ICC, ELISA
Cross Reactivity	Human, Mouse, Rat
SWISS	Q13490
Protein Weight	70kDa
Shipping	Ice bag