


Anti-Ciapiin1 Polyclonal Antibody For Wb, If/Icc - Q6Fi81

Item Number: CM0029187



Introduction

CIAPIN1 rabbit polyclonal antibody for WB, IF/ICC, IP, ELISA applications. Recognizes human and mouse CIAPIN1, a Fe-S cluster assembly protein involved in cytosolic iron-sulfur biogenesis and apoptosis regulation. Optimal for research in cell death and iron metabolism.

[Learn More](#)

Parameter	Value
Product Name	CIAPIN1 Rabbit pAb
Remarks/Alias	DRE2; CIAE2; PRO0915; Anamorsin; CIAPIN1
Species	Human
GeneID (Human)	57019
GeneID	57019
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1-312 of human CIAPIN1 (NP_064709.2)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IF/ICC, IP, ELISA
Cross Reactivity	Human, Mouse
SWISS	Q6Fi81
Protein Weight	34kDa
Shipping	Ice bag