

Anti-Cic (Capicua) Rabbit Polyclonal Antibody For Wb And Elisa - Q96Rk0

Item Number: CM0013996



Introduction

Rabbit polyclonal antibody targeting human CIC (capicua homolog), suitable for Western blot and ELISA applications. Cross-reacts with mouse. Useful for studies of CNS development and transcriptional repression.

[Learn More](#)

Parameter	Value
Product Name	CIC Rabbit pAb
Remarks/Alias	MRD45; CIC
Species	Human
Gene ID (Human)	23152
Gene ID	23152
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 1379-1608 of human CIC (NP_055940.3)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q96RK0
Protein Weight	164kDa
Shipping	Ice bag