


Anti-Cidea Rabbit Polyclonal Antibody For Wb, Ihc-P, Elisa - O60543

Item Number: CM0025371



Introduction

Anti-CIDEA rabbit polyclonal antibody validated for WB, IHC-P, and ELISA. Reacts with human, mouse, and rat. Targets lipid droplet fusion and apoptosis. Ideal for metabolic research.

[Learn More](#)

Parameter	Value
Product Name	CIDEA Rabbit pAb
Remarks/Alias	CIDE-A; CIDEA
Species	Human
Gene ID (Human)	1149
Gene ID	1149
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 20-219 of human CIDEA (NP_001270.1)
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, IHC-P, ELISA
Cross Reactivity	Mouse, Rat
SWISS	O60543
Protein Weight	25kDa
Shipping	Ice bag