


Anti-Cideb Polyclonal Antibody For Wb, Elisa - Q9Uhd4

Item Number: CM0019387



Introduction

Rabbit polyclonal antibody against human CIDE-B, validated for WB and ELISA. Detects endogenous CIDE-B in lipid droplet studies, cross-reacts with mouse. Suitable for hepatocyte biology and VLDL maturation research.

[Learn More](#)

Parameter	Value
Product Name	CIDE-B Rabbit pAb
Remarks/Alias	CIDE-B; cell death activator CIDE-B; IN; LHR; MC56; MDU2; MDU3; MIC4; Pgp1; CDW44; CSPG8; HCELL; HUTCH-I; ECMR-III
Species	Human
GeneID (Human)	27141
GeneID	27141
Immunogen	Recombinant fusion protein containing a sequence corresponding to amino acids 80-160 of human CIDE-B (NP_001305736).
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Mouse
SWISS	Q9UHD4
Protein Weight	25kDa
Shipping	Ice bag