


Anti-Cideb Rabbit Polyclonal Antibody For Wb, Elisa - Q9Uhd4

Item Number: CM0021272



Introduction

Rabbit polyclonal antibody against human CIDE-B, validated for WB and ELISA. CIDE-B is a lipid transferase involved in lipid droplet fusion, VLDL maturation, and HCV entry. Suitable for liver and metabolic disease research.

[Learn More](#)

Parameter	Value
Product Name	CIDE-B Rabbit pAb
Remarks/Alias	CIDE-B; cell death activator CIDE-B; IN; LHR; MC56; MDU2; MDU3; MIC4; Pgp1; CDW44; CSPG8; HCELL; HUTCH-I; ECMR-III
Species	Human
GeneID (Human)	27141
GeneID	27141
Immunogen	Recombinant Protein
Source	Rabbit
Category	Polyclonal Antibodies
Application	WB, ELISA
Cross Reactivity	Human
SWISS	Q9UHD4
Protein Weight	25kDa
Shipping	Ice bag